



---

# ***Agility.....***

---

KYMCO Healthcare's new ***Agility*** 8mph -12.8 km/h mid-size mobility scooter has been designed to provide you with the ultimate in performance, luxury, comfort & manoeuvrability. Designed using the very latest technologies, ***Agility*** will delight & surpass your mobility expectations!

With stylish large 12-inch -305 mm low profile alloy wheels & a full oil damped suspension system, this new model will provide you with the ultimate in comfort. Even when driving in the more challenging tough environments you will enjoy a trouble free & comfortable journey.

We have included many new high-end features for your convenience on the new ***Agility***, as standard !

## **Features Include :**

- Waterproof backlit LCD display – Leading motorcycle technology.
- Easy to use delta handlebar with high visibility rear view mirrors.
- Soft touch illuminated waterproof switches for safe easy use.
- Luxury breathable fully adjustable Captain Seat for extra comfort.
- Convenient 12v power socket – You can charge your mobile phone onboard.
- Low consumption LED front & rear lighting system.
- Easy access on board storage facility for convenience & safe storage.
- The latest modern design lines with 5 attractive metallic colours to choose from.
- Automotive style dashboard with a service interval warning.
- Superior ground clearance for outdoor driving.
- Compact & easy manoeuvrability means you can drive where you want to go.

***With Kymco's superior reliability, quality & design the brand new Agility is the future of mid - size luxury mobility travel !***

***Why not try one out today & visit your local  
KYMCO Healthcare Retailer !***

Quality, Reliable, Stylish Luxury !



Waterproof LCD Dash



Oil Suspension



LED Lighting



12" Low Profile Alloys

### Agility Technical Data

<b>Spfications</b>	<b>EQ 35FA (UK)</b>
Maximum User Weight .....	165 kg - 26 Stone
Overall Length.....	1,280 mm
Overall Width .....	575 mm
Height .....	1,100 mm
Max Speed .....	12.8 km/h - 8 mph
Travel Range .....	40 km - 25 miles
Wheel Size .....	12" - 305 mm
Wheel Type .....	Low Profile Pneumatic Alloy
Ground Clearance .....	138 mm - 5.4"
Total Weight (inc batteries) .....	107 kg
Batteries .....	50 ah
Seat Height - From Deck .....	460 mm - 520 mm
Lights Fitted .....	Yes (LED)
Turning Radius .....	1,560 mm
Maximum Slope .....	12 deg

- Glossy Bronze
- Metallic Mink
- Cherry Red
- Glossy Black
- Sapphire Blue



DEALER STAMP

**KYMCO Healthcare UK Limited**  
**30 North Road**  
**Bridgend Industrial Estate**  
**Bridgend CF31 3TP**

**t: +44(0) 1656 645755 f: +44(0) 1656 661983**

**e: [sales@kymcohealthcare.co.uk](mailto:sales@kymcohealthcare.co.uk)**

**w: [www.kymcohealthcare.com](http://www.kymcohealthcare.com)**

 **KYMCO Healthcare UK**  
**@KYMCO\_ForU**



Some product images may include non-standard accessories for illustrative purposes.  
Product details are correct at time of printing and may change without prior notice.