

2015/16

ottobock.

Item Price - £1258.00	Order no.: -----
-----------------------	---------------------

Kimba Cross Paediatric postural positioning system

Maximum load capacity of the outdoor mobility base: 50 kg/110 lbs (including seating unit)
Article number: HR3401=0000_K

Maximum load capacity of the seating unit: 40 kg
Article number: 470G71=20000_K

<input type="checkbox"/> Provide estimate only	<input type="checkbox"/> Order
<input type="checkbox"/> Attachment for custom fabrication enclosed	
Billing address/customer number	
Delivery address	
Sales person	Name Date Signature



Standard equipment

- Three-wheeled outdoor mobility base with telescoping frame, for mounting a Kimba neo seating unit, size 1 or size 2, or a custom seating shell
- Can be folded/disassembled
- Elastomer spring system
- Drum brake
- 18" rear wheels with pneumatic tyres
- Rigidly mounted 16" front wheel with pneumatic tyre, includes splash guard
- Height-adjustable push bar

For children from approx. 1 to 10 years

Outdoor mobility base		
	<input type="checkbox"/>	KIMBA Cross outdoor mobility base £1258.00
Frame colour		
016	<input type="checkbox"/>	Silver metallic
013	<input type="checkbox"/>	Jet black RAL 9005 matte
For use in conjunction with		
	<input type="checkbox"/>	Kimba Neo seating unit, size 1
	<input type="checkbox"/>	Kimba Neo seating unit, size 2
	<input type="checkbox"/>	Mygo Seat seating unit
	<input type="checkbox"/>	Squiggles Seat seating unit
	<input type="checkbox"/>	Lars car seat, size 1
	<input type="checkbox"/>	Lars car seat, size 2

Seat adapter for Kimba neo seating unit		
HR34011300	<input type="checkbox"/>	Size 1
HR34021300	<input type="checkbox"/>	Size 2
HR32x48300	<input type="checkbox"/>	Seating shell interface

Accessories for outdoor mobility base		
HR34012100	<input type="checkbox"/>	12" swivelling wheel with quick-lock attachment (will be delivered in addition to 16" front wheel) £164.00
HR34017000	<input type="checkbox"/>	Anti-tipper, complete £114.00
HR34018000	<input type="checkbox"/>	Clothing protector (plug-on wheel covers) £164.00
HR34018100	<input type="checkbox"/>	Clothing protector extension around the wheel (not available in combination with spoke protector) £62.00
HR34019000	<input type="checkbox"/>	Spoke protector, transparent £59.00
HR34015000	<input type="checkbox"/>	Storage bag (e.g. for the swivelling front wheel) £49.00
HR34016100	<input type="checkbox"/>	Bicycle draw-bar, complete (includes coupling from Weber company, type E, safety strap)* £205.00
HR34016200	<input type="checkbox"/>	Equipment package required according to German Road Traffic Licensing Regulations (StVZO) (lighting/pennant/roll bar) £106.00
HR16001102	<input type="checkbox"/>	"Techno" buggy board £101.00
HR32xx2700	<input type="checkbox"/>	Footrest assembly for seating shell interface £270.00

* According to a leaflet for bicycle trailers issued by the German Federal Ministry of Transport: (www.pdeleuw.de/fahrrad/fakten) – Recommended overall weight of the bicycle trailer incl. load: 40 kg. Recommended in combination with four-point belt, groin strap, spoke protectors.

470G71=20000_K Kimba Neo seating unit		
KC 01	<input type="checkbox"/>	Size 1 (from approx. 1 to 6 years) £838.00
KC 02	<input type="checkbox"/>	Size 2 (from approx. 4 to 10 years) £919.00

Frame colour		
-016	<input type="checkbox"/>	Silver metallic Inclusive
-013	<input type="checkbox"/>	Jet black RAL 9005 matte Inclusive

Padding		
KD 00	<input type="checkbox"/>	Without Inclusive
905	<input type="checkbox"/>	Spacer fabric beige/brown ¹⁾ Inclusive
906	<input type="checkbox"/>	Spacer fabric beige/green ¹⁾ Inclusive
KD 09	<input type="checkbox"/>	Microfibre grey/berry £111.00
KD 21	<input type="checkbox"/>	Microfibre grey/blue £111.00
KD 22	<input type="checkbox"/>	Microfibre black/berry £111.00
KD 23	<input type="checkbox"/>	Microfibre black/blue £111.00

Back insert		
KD 01	<input type="checkbox"/>	Standard Inclusive
KD 02	<input type="checkbox"/>	With integrated headrest £185.00

Footrest		
KB 03	<input type="checkbox"/>	Standard Inclusive
KB 05	<input type="checkbox"/>	Angle-adjustable £76.00

Footrest accessories			
KB 08	<input type="checkbox"/>	Anti-slip mat	£19.00
KB 09	<input type="checkbox"/>	Padding	£52.00
KB 06	<input type="checkbox"/>	Foot straps, pair	£21.00
KB 07	<input type="checkbox"/>	Hook-and-loop foot straps with instep pad	£24.00
Seating unit accessories			
KD 03	<input type="checkbox"/>	Headrest supports ¹⁾ (in padding colour)	£148.00
KD 04	<input type="checkbox"/>	Thoracic supports ¹⁾ (in padding colour)	£148.00
KP 02	<input type="checkbox"/>	Lap belt, padded ¹⁾	£55.00
KP 03	<input type="checkbox"/>	Four-point belt, padded ¹⁾	£85.00
KP 05	<input type="checkbox"/>	Five-point belt with small buckle (recommended for seat, size 1)	£59.00
KP 04	<input type="checkbox"/>	Five-point belt, large buckle	£59.00
KP 06	<input type="checkbox"/>	Groin strap	£91.00
KP 07	<input type="checkbox"/>	Pelvic band, size 1 (only for seat, size 1)	£107.00
KP 08	<input type="checkbox"/>	Pelvic band, size 2	£107.00
KP 09	<input type="checkbox"/>	Pelvic band, size 3 (only for seat, size 2)	£107.00
KP 10	<input type="checkbox"/>	Fixation vest	£107.00
KP 11	<input type="checkbox"/>	Chest and shoulder support, size 1	£107.00
KP 12	<input type="checkbox"/>	Chest and shoulder support, size 2	£107.00
KP 13	<input type="checkbox"/>	Chest and shoulder support, size 3	£107.00
KQ 01	<input type="checkbox"/>	Grab rail	£47.00
KQ 02	<input type="checkbox"/>	Padding for grab rail (in padding colour)	£30.00
KQ 03	<input type="checkbox"/>	Grab rail with PVC cover	£54.00
KE 01	<input type="checkbox"/>	Armrests	£79.00
KC 06	<input type="checkbox"/>	Abductor, size 1 (L 9 cm x W 6 cm x H 9 cm), only for seat, size 1, in padding colour	£99.00
KC 07	<input type="checkbox"/>	Abductor, size 2 (L 10 cm x W 7.5 cm x H 10 cm), in padding colour	£99.00
KC 08	<input type="checkbox"/>	Abductor, size 3 (L 12 cm x W 9 cm x H 10 cm), only for seat, size 2, in padding colour	£99.00
KD 10	<input type="checkbox"/>	Lumbar support padding	£30.00
KD 11	<input type="checkbox"/>	Neck support	£43.00
KS 01	<input type="checkbox"/>	Tray, grey	£179.00
KS 02	<input type="checkbox"/>	Tray, transparent	£253.00
KD 12	<input type="checkbox"/>	Incontinence cover	£106.00
KW 01	<input type="checkbox"/>	Summer-winter slipsack, size S	£231.00
KW 02	<input type="checkbox"/>	Summer-winter slipsack, size M	£231.00
KW 03	<input type="checkbox"/>	Summer-winter slipsack, size L	£231.00
KS 08	<input type="checkbox"/>	Backpack, attaches to push bar	£36.00
KS 04	<input type="checkbox"/>	Baby Changing bag, brown*	£101.00
KS 05	<input type="checkbox"/>	Baby Changing bag, green*	£101.00
KS 06	<input type="checkbox"/>	Baby Changing bag, black*	£101.00
KS 07	<input type="checkbox"/>	Cupholder	£19.00
KW 04	<input type="checkbox"/>	Rain cape, grey	£53.00

KW 05	<input type="checkbox"/>	Canopy in padding colour, incl. adapter and rain cover ¹⁾	£195.00
KW 06	<input type="checkbox"/>	Rain canopy, transparent, with cover for the legs, includes adapter	£137.00
KP 01	<input type="checkbox"/>	Selected accessories mounted (only applies to accessories marked with ¹⁾)	£65.00
KC 10	<input type="checkbox"/>	Seat depth shortener	£67.00
KD 20	<input type="checkbox"/>	Back angle limitation	£69.00
KP 14	<input type="checkbox"/>	Mounting of installation kit (affixing clamps to the back plate etc.)	£71.00
*Contains a vast range of accessories such as: changing mat, soothers bag, insulated bottle bag and a plastic utensil bag, as well as numerous compartments and zippered pockets.			

Notice: Seating unit, size 2, is available in custom fabrication upon request.

HR34016200 – Equipment package required according to German Road Traffic Licensing Regulations (StVZO)

According to a leaflet for bicycle trailers issued by the German Ministry of Transport: (www.pdeleuw.de/fahrrad/fakten)

Recommended overall weight of the bicycle trailer, including load: 40 kg. Recommended in combination with four-point belt, groin strap, spoke protectors. Seating shell interface and footrest assembly in combination with Squiggles and Mygo not possible.

HR34018000 – Clothing protector is a plug-on wheel cover

HR34018100 – Clothing protector extension around the wheel cannot be used in combination with spoke protector

Accessories only for corresponding size

Footrest assembly only available in combination with seating shell interface