

Otto Bock®

QUALITY FOR LIFE



Eco-Buggy

Your Child Deserves The Best

Sturdy and Stable – Aluminium frame is light yet rigid for comfort, positioning and safety. Up to 50 kg (110 lbs) User weight. Reinforced nylon seat back fabric provides support and comfort.

Lightweight – Only 9 kgs (20 lbs) Compact Folds to a compact size for easy manoeuvrability and storage. Handy storage pocket under seat

Comfort and Flexibility – Adjustable backrest (90° to 110 °) and foot-rest (20 - 34 cm) allowing custom comfort settings with adjustment for your child's growth.

Versatile Swivel front wheels with locking mechanism for easier transit over rough terrain. Sprung rear wheels with brake ensures user safety and comfort.

Innovation and Quality Traditional Otto Bock design and build craftsmanship together with modern design and manufacturing assures you of the very best for your child.



Otto Bock HealthCare plc

32, Parsonage Road · Englefield Green GB-Egham, Surrey TW20 0LD · Phone +44 1784-744900 · Fax +44 1784-744901
bockuk@ottobock.com · www.ottobock.co.uk



Eco-Buggy



Sturdy seat and back upholstery



Double castors with front steering restraint



Reclining back positioning



Back wheels with locking mechanism



Height adjustable footrest



Folding storage pocket

Technical Data

Seat Depth:	30 cm
Seat Width:	35 cm
Back Rest Height:	55 cm
Back Rest Angle:	90° - 110°
Distance from seat edge to footrest	20 - 34 cm*
Weight:	9 kg
Weight Capacity:	50 kg
Measurement when folded:	30 x 40 x 12
Wheel Ø front/back:	16.5 / 20

* max. lower leg length requires adjustment of the casters



Frame Colour



anthracite metallic

Nylon Upholstery



blue

Otto Bock®

QUALITY FOR LIFE