



K-ACTIV

The new **KYMCO K-ACTIV** Powered Wheelchair was designed to be the most flexible fixed frame powerchair, providing you with the ultimate in comfort, performance & independence! With multiple adjustment parameters this compact chair can be adjusted to suit your specific needs. **K-ACTIV** is also extremely maneuverable allowing you to drive with confidence both outdoors & indoors with ease. You can go where you want with **K-ACTIV**!

Some of the key features that will benefit you are listed below:

- Adjustable rear suspension for a smooth & comfortable drive
- Fully programmable control system - Set the performance to suit how you want to drive
- 14"- 35cm Pneumatic rear drive wheels for smooth comfortable driving
- 9"- 23cm Puncture proof front tyres for peace of mind on your journey
- Comfort seat cushion & tension adjustable backrest upholstery for extra comfort
- Easy to adjust, angle adjustable backrest for a perfect driving position
- Multi-adjustable seating system with width & depth adjustments
- Independently adjustable joystick to suit your driving position
- 6 mph / 10 kph class 3 - sporty driving for outdoors use
- 50 ah batteries - 27mile / 44km, extra long driving range
- Powerful 450w motors for driving in tough environments
- Low consumption LED lighting system
- Available in attractive Metallic Mink & Golden Lemon



Golden Lemon



Metallic Mink

Specifications:

Overall Length	108 cm
Overall Width	68 cm
Overall Height	105.5 cm
Unit Height (backrest folded down)	62.5 cm
Total Weight (inc'g Batts)	97 kg
Total Weight (excl'g Batts)	66 Kg
Wheels Size Rear	35 cm / 14"
Wheel Size Front	23 cm / 9"
Seat Width (adj arms)	46 - 64 cm (18" - 25")
Seat Depth (adj back/footrest brackets)	46 - 56 cm (18" - 22")
Adjustable seat tilt angle	0, 2, 5, 8, 10 deg
Backrest Height	51 cm (20")
Seat to Floor Height	49.5 cm (19.5") std 5 deg tilt
Footrest	Swing in / Swing out
Footplate	Angle Adjustable & heel loops
Turning radius	55.5 cm (22")
Ground Clearance	12.5 cm (5")
Ground Clearance of anti tipper wheels	6.5 cm (2.5")
Max Speed	6 mph / 10 km/h
Max User weight	25 stone / 158kg
Max Travel Range	27 miles / 44 Km
Battery Size (12v x 2)	50 Ah
Manual Handbrake	Yes
Suspension	Adjustable
Charger	24v 5A
Motor Size	450 W
Max Safe Slope	10 deg



**CRASH TESTED
ISO 7176-19**



KYMCO Healthcare UK Limited
30 North Road,
Bridgend Industrial Estate,
Bridgend CF31 3TP

t: +44(0) 1656 670095 **f:** +44 (0) 1656 858353

e: sales@kymcohealthcare.co.uk

w: www.kymcohealthcare.com

DEALER STAMP




Swingaway Joystick



Adjustable footrest



Led Lighting



Rear Suspension



Product details are correct at time of printing and may change without prior notice!