

Aquatec® | **Orca Bath Lift Range**



! Hand Control
Battery contained
in the hand set



3 YEAR WARRANTY



Aquatec **Bath Lifters - Tried and trusted by healthcare professionals.**

The *Aquatec Orca* family is comprised of superior quality and durable bath lifters that are designed to offer comfort, safety and reassurance to those most in need of support as they bathe. The *Orca* incorporates increased leg room, ease of disassembly and an ergonomic, floating hand control that contains the light weight battery. The *Orca* has an impressive weight limit of 135 kg (21.5 stone) as standard, while the *Orca XL* will support patients up to 170 kg (26.5 stone).

An emphasis on ergonomic design allows the *Orca* to sit lower and further towards the back of the bath, giving increased leg room for deeper and more relaxed bathing. The increased height of the backrest provides full support whilst upright or reclined.

AQUATEC®



Yes, you can.®



1470754 : Aquatec Orca

with reclining backrest

- Quick release hand control
- Battery with Sony™ technology
- Auto unlock for quick removal from bath
- **1470754 - White**
- **1470764 - Blue**



Aquatec Orca F

An upright bath lift with a fixed backrest, that offers support for those who prefer to bathe upright. The integrated head support of the high backrest ensures a safe and comfortable bath.

- **Aquatec Orca F** with fixed backrest
- **1471074 - White**



Aquatec Orca XL

A powerful bath lifter with maximum client weight of 170 kg (26.5 stone).

The reinforced scissor mechanism combines safe lifting with a compact construction. Minimum seat height of just 60 mm; allows plenty of room for a deeper and more comfortable bath.

- **Aquatec Orca XL** with reclining backrest
- **1471068 - White**

Features and Options



Compact hand control with integral battery

Orca bath lifts incorporate a revolutionary lightweight battery contained within the ergonomic hand control. Raised profile, coloured buttons and an emergency stop feature allow simple, controlled operation to give reassurance and independence.



Safe hand control fitting

The hand control connects within seconds to the back of the Orca, providing a secure and completely watertight seal. Effortlessly released with just one hand.



Two-piece design

When pushed forward the **Orca's** backrest folds down and automatically unlocks. The backrest and base section can then be removed from the bath with ease. The 2 sections fit neatly into one another to give an extremely compact size for storage.



Textured seating surface

The **Orca** has a textured surface allowing use without covers if desired. The side flaps are hinged flush to the main seating surface and allow comfortable and safe side transfers, without entrapment of skin.



Supported transfers

For those who require additional support during transfer, the **Aquatec Trans** allows the client to slide onto the seat and swivel into the bath with the minimum of effort.



Refresh kit

The **Orca** refresh kit comprises a full set of covers, oyster suction cups, side flap diverters and user guide. Ideal for re-issue



Height adaptors

The **Orca** can be adapted to suit a full range of bath heights. Height adaptors can be quickly fitted to the bath lift base.



Large oyster cups

These extra large oyster cups allow the **Orca** to adhere to steeply curved and textured baths providing exceptional stability and safety. Smaller cups are available as an option.

Features and Options



Transport bag

A sturdy canvas bag with wheels for discreet transport and safe storage.



Non-slip cover mats

Cover mats are available in white or blue and are easily removed for machine washing at 60°C.



Fast Fit Pommel

A cost effective and easy to fit option (available in 7cm and 12cm diameter in blue colour only).



Easy to clean

The open, clean-cut design of the **Orca** combines with its smooth surfaces to offer effortless cleaning.

Technical data



Orca/Orca F	36 cm	135 kg
Orca XL	36 cm	170 kg

Colours



White



Blue

INVACARE® Limited
 Pencoed Technology Park - Pencoed -
 Bridgend CF35 5HZ
 Tel: +44 1 656 776 222
 Fax: +44 1656 776 220
 E-mail: uk@invacare.com
www.invacare.co.uk
 Sales Order E-mail:
ordersuk@invacare.com

INVACARE® Ireland Ltd
 Unit 5 - Seatown Business Campus
 Seatown Road - SWORDS - County
 Dublin
 Tel. +353 1 8107084
 Fax +353 1 8107085
 Email: Ireland@invacare.com
www.invacare.ie

© Invacare® International Sarl 2010
 All rights reserved.

All measurement have been done under the Invacare® standard test conditions, unless quoted differently. All information quoted is believed to be correct at time of print.
 Invacare® reserves the right to alter product specifications without prior consultation.

Orca Family - UK - 09/2010 - 1503848



Yes, you can.®