

Pride QUEST



A Light, Portable and
Easy to Fold Scooter



PRIDE® MOBILITY SCOOTERS
LIVE YOUR BEST®

Pride®
Mobility Products Ltd.
— LIVE YOUR BEST® —



Tel: 01869 324600
Email: sales@pride-mobility.co.uk
Web: www.pride-mobility.co.uk

Features and Benefits

Introduction

The Pride® Quest Folding Scooter is perfect for the active lifestyle. Offering great performance indoors and out at a maximum speed of 4mph. The Quest provides a comfortable ride and a simple easy-fold system. The Quest is portable and adjusts to fit perfectly inside any small space. The Quest is sleek, stylish and extremely easy to operate.

Front and Rear Lights

Featuring highly visible front and rear LED lights, which have been expertly incorporated within the frame of the Pride® Quest.

These LED lights can be accessed from the delta tiller console. Once switched on they enable the Quest to be visible to on-coming vehicles and public.



Flexible Quick



Front LED Lights



Rear LED Lights

Features and Specification

Detachable Battery Box

Each Pride® Quest comes with a single 11.5 AH Lithium-Ion battery pack. This pack can be easily detached from the scooter, and features an anti-theft key lock.



Available Colours

The Pride® Quest comes in two optional colours:














Red



Black

Folding Scooter

Technical Specifications

Weight Capacity	Max Speed	Ground Clearance	Turning Radius	Overall Length	Overall Width	Front/Rear Tyres	Total Weight w/o batteries	Heaviest piece	Standard Seating	Batteries Required
 18 Stone (115 kg)	 4 Mph (6.4 km/h)	 6 cm (2.36 inch's)	 82 cm (32.2 inch's)	 93 cm (36.6 inch's)	 48.5 cm (19.09 inch's)	 Front: 20cm Rear: 23 cm	 3.6 Stone 23 k.g.	 N/A	 Folded Seat 12.6 k.g.	 x1 11.5 Ah LI-ION battery

Warranty

Two-year limited warranty on frame
One-and-a-half-year limited warranty on Controllers
One-year limited warranty Electronic Components, battery and Charger.

Notes

Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure.
Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.
Battery weight may vary based on manufacturer.

Tel: 01869 324600
Email: sales@pride-mobility.co.uk
Web: www.pride-mobility.co.uk

Pride[®]
Mobility Products Ltd.
— LIVE YOUR BEST[®] —

Pride[®]

Mobility Products Ltd.

— LIVE YOUR BEST —



Please scan
code for
more information

*images used for illustration purposes only
(Version 1: April 2015)

Tel: 01869 324600
Email: sales@pride-mobility.co.uk
Web: www.pride-mobility.co.uk