



# Tilt-in-Space Rodeo



















## Standard Specifications

### MEASUREMENTS:

INCHES

HCPCS:E1232

\* - Plus or minus .5 inches due to fabric  
 Δ - Weight of chair is without front rigging

	RD 10	RD 12	RD 14	RD 16
Seat Width* 	10	12	14	16
Seat Depth* 	7 – 12	9 – 14	11 – 16	13 – 18
Seat Back Height 	21	22	24	25
Seat to Footplate* 	9 – 14	9 – 14	11 – 16	12 – 18
Seat to Floor 	20	21	22	22
Recline Adjustment 	0°, 20°	0°, 20°	0°, 20°	0°, 20°
Tilt Adjustment 	5° - 45°	5° - 45°	5° - 45°	5° - 45°
Overall Height 	43	43	44	45
Overall Length 	30.5	35	35	36
Overall Width 	21	23	25	27
Folded Length 	33.5	34	34	40
Folded Width 	21	23	25	27
Folded Height 	17.5	17	18	18
Headrest Extension 	6	8	8	8
Shoulder Strap Height 	13 – 19	15 – 21	16 – 22	18 – 24
Footplate Adjustment Angle 	±20°	±20°	±20°	±20°
Weight Capacity Regular/Transit 	75/66 lbs	75/66 lbs	100/100 lbs	170/140 lbs
Weight of Chair Δ 	29 lbs	33 lbs	37 lbs	39 lbs





















# Tilt-in-Space Rodeo

## Standard Specifications

### MEASUREMENTS:

CENTIMETERS

\* - Plus or minus 1.25 cm due to fabric  
 Δ - Weight of chair is without front rigging

	RD 10	RD 12	RD 14	RD 16
Seat Width* 	25.5	30.5	35.5	40.5
Seat Depth* 	17.8 – 30.5	22.9 - 35.6	27.9 - 40.6	33 - 45.7
Seat Back Height 	53.3	56	61	63.5
Seat to Footplate* 	23 – 35.6	23 – 35.6	28 - 40.6	30.5 – 45.7
Seat to Floor 	51	53	56	56
Recline Adjustment 	0°, 20°	0°, 20°	0°, 20°	0°, 20°
Tilt Adjustment 	5° - 45°	5° - 45°	5° - 45°	5° - 45°
Overall Height 	109	109	112	114
Overall Length 	77.5	89	89	91.5
Overall Width 	53	58.5	63.5	68.5
Folded Length 	85	86	86	101.5
Folded Width 	53.3	58.5	63.5	68.5
Folded Height 	44.5	43	45.7	45.7
Headrest Extension 	15	20	20	20
Shoulder Strap Height 	33 - 48.3	38 - 53.3	40.6 - 56	45.7 – 61
Footplate Adjustment Angle 	±20°	±20°	±20°	±20°
Weight Capacity Regular/Transit 	34/30 kg	34/30 kg	45.5/45.5 kg	77/63.6 kg
Weight of Chair Δ 	13.2 kg	15 kg	16.8 kg	17.7 kg